



HWS Fermenter
Economy 20:

Modular construction
guarantees the individual
fulfillment of your process
conditions.

H_oW_oS_o Economy Fermenter 10/20

Contents

1) Characteristics

2) Construction and Function

3) Control System

3.1 General

3.2 Thermostatisation

3.3 Stirring

3.4 pH Measurement /- control

3.5 pO₂ Measurements

4.) Technical Specifications

5.) References

Short Characteristics

The HWS Economy fermenter is giving you the highest flexibility. You can desire:

different working volumes 1 , 2 , 2.5 , 5 , 6 , 10 and 20 ltr

stirrer systems, motor with/without PC-connection, different stirrer executions

temperature variation by circulators

measuring and control by different modules

You can adapt the HWS Economy fermenter to specific process conditions to have highest employment.

The fermenter can be sterilized in an autoclave.

The system is easy to handle, all parts can be used even by un-trained personal.

The control units are equipped with outputs for registration of data.



Construction and Function:

The HWS- Economy Fermenter is available in 2 different versions:

Economy 10 fermenters with working capacities from 1-5 ltr are executed as desk systems

Economy 20 fermenters with capacities from 6-20 ltr are installed in a moveable stainless steel frame

The culture vessels are manufactured from borosilicate glass tubes with strong walls and a thermo jacket, for good temperature controlling.

The flanges contain a groove for a simple insert of O-ring seals.

The covers are manufactured either in glass (Economy10) or in stainless steel (Economy20).

The stirrers are sealed

Economy 10- dynamic seal

Economy 20- mechanical seal

All covers are equipped with sufficient ports for sensors feeding and harvesting:

Economy 20- 12 ports (stainless steel)

Economy 10- 07 ports (glass,plastic)



The ports are executed in screw threads and caps with septa, or closing caps with seals. In use this has shown to be most practicable.

The stainless steel stirrer shaft can be equipped with different stirrer systems like disc blade, propeller and others.

The systems are easy to exchange and fixable.

The temperature regulation of the culture vessel is achieved by a circulator, which guarantees a high temperature accuracy. Therefore it can be used for thermophile organisms.

The supply of air in the fermenter can be produced by compressed cleaned air with a valve and a gas flow meter; while the air introduction is efficiently aimed by a PTFE distribution ring in the medium.

Other parts like control systems and peristaltic pumps are completing the systems.

Control System

3.1 General

The measuring and control of the process can be executed by the use of different modules.

These modules are mounted in a 19" rack, containing the main power supply.

The different modules:

- pH measurement/control
- pO₂ measurement
- pO₂ control
- dosing pumps for acid/alkaline
- valve for CO₂
- PT 100 for temperature registration

Stirrer motor and circulator can be equipped with outputs for giving their data to a suited printer.

3.2 Thermostatisation

The heating of media liquids results by a closed circulator and the thermo jacket of the glass vessel. The temperature can be measured by an external PT 100.

The prescribed temperature can be inserted by an adjustment.



3.3 Stirring

The stirrer system contains a strong motor which enables a constant rotation. The speed is indicated by a digital display. If desired the stirrer motor can indicate torque. Speed and torque can be registered by a 0-20 mA output.



Control System

3.4 pH Measurement/-control

The measurement of the pH-values is executed by a sterilizable electrode by Schott. The cable is connected to the HWS pHR 1402 device, which has to be calibrated before use. The accuracy enables the connection to dosing pumps. By this constellation a regulation of the pH value is guaranteed.



3.5 pO₂ Measurement

The measurement of the solved oxygen is achieved by a oxygen transmitter 4050e. Each transmitter is equipped with a micorprocessor offering powerfull features through a user-friendly menu-driven Setup programm.

With a potentiometer the zero value can be compensated and with a second measurement in oxygen saturated solution the gradient can be fixed.

The instrument also supports automatic temperature compensation from 0...50°C.

Measurement in ppm, mg/l or % saturation.

All measured data can be shown later with a printer.



Technical Specifications

Culture vessel

Working Volumes	ltr.	1 2 2.5 5	6 10	20
Diameter int.	mm	150	200	260
Cover		Borosil.glass	s/s	s/s
Ports		7	12	12
Seals		viton/silicone	viton/silicone	viton/silicone

Circulator:

Sensor	PT 100
Working range	0...200°C
Accuracy	0.1°
Constant	0.05°
Heating	2kW
Cooling	water
Control	external PT100
Power supply	230V/50Hz

Stirrer motor:

Power output	30-100W
Speed	50...2000 rpm
Indicator digital display LCD	
in special execution: with torque indicator and values output	
	4-20mA
Power supply	230V/50 Hz

Technical Specifications

pH measurement HWS pHR1402:

Sensor:	Sterilizable electrode from Schott
Range:	1-13 pH
Temperature compensation:	automatically
Output:	4-20 mA

Control:

Range:	1 - 13
Adjustability:	0.01
Constancy:	0.01

Temperature:

Temperature indicator:	integrated in HWS pHR 1402
Resolution:	0.01°
Output:	4-20mA

pO₂ measurement of solved oxygen:

Resolution	0.01 mg/l or 0.1 %
Accuracy	+/- 1% of full scalereading
Sensor po ₂	polarographic type (Clark)
Sensor Temp.	NTC 22k ohm thermistor
Temperature compensation	auto/manual
Output	4...20mA
Power supply	80...250 VAC

